

Qosina Large Bore Touhy Borst Adapter with Side Port

Edgewood, NY, January 24, 2005 – Qosina is pleased to announce a new custom-designed addition to its kit components, the Large Bore Touhy Borst adapter with side port and threaded flare connector.

- ◆ **Part #11182** Large Bore Touhy Borst Adapter with threaded flare connector, red, 7-17 FR
- ◆ **Part #11183** Large Bore Touhy Borst Adapter with threaded flare connector, blue, 0-13 FR

This adapter is available in two sizes, 0-13 French and 7-17 French. This part was designed to Qosina's requirements as an addition to our existing Touhy Borst (Part #s 22827, #22828) adapter selection. This side port accepts 5mm OD or 7.7mm ID size tubing to allow an easy access to the fluid path.

This new adapter has a custom molded tapered, silicone gasket to provide a tight seal with the adjustable, easy-to-control cap. The shape, size and style are enhanced for additional versatility. The polycarbonate body provides great clarity.

For additional kit components, please visit our comprehensive website, qosina.com or call our Customer Specialists for help with your sampling/ catalog request or order.

Qosina is an ISO 9001:2000 registered worldwide supplier of stock and custom-made OEM disposable components including swabs and applicators, stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing. Custom sourcing service is available by knowledgeable specialists using Qosina's extensive resources.

