



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

QOSINA STRAIGHT CONNECTOR

Edgewood, NY, September 21, 2004 – Qosina offers economically priced straight connectors used primarily in the construction of tracheal tube sets. These straight connectors (**#31300 to #31316**) are sized from 15 mm x 2.5 mm, in 0.5 mm increments, to 15 mm x 11.5 mm. The straight connectors have all polysulfone construction for clarity.

Qosina offers additional bulk, non-sterile tracheal kit components including tracheal airway connectors (#53519), tracheal circuit connectors (#53525) and double swivel elbows (#56025).

For additional information regarding these products, please visit our comprehensive website, qosina.com or call our Customer Specialists for help with your sampling/ catalog request or order.

Qosina is an ISO 9001:2000 registered worldwide supplier of stock and custom-made OEM disposable components including swabs and applicators, stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, beakers, gloves, syringes and tubing. Custom sourcing service is available by knowledgeable specialists using Qosina's extensive resources.

