



NEWS RELEASE  
FOR IMMEDIATE RELEASE

## Qosina Extra Large Bore High Flow Threaded Connectors

Edgewood, NY, February 22, 2005 – Qosina is pleased to present our **extra large bore high flow threaded connectors**. These high quality connectors allow for faster delivery of viscous fluids and avoid the risk of connection errors between incompatible fluid systems. The extra large bore threaded male connector (**Part # 98663**), and female connector (**Part # 98633**) are composed of ABS, one and one half inches in length and have a 6mm OD.

The 5mm large, 6mm extra large and 8mm super large bore threaded connectors will only mate up to each other and not to standard luers. Other Qosina components that can be used in conjunction with our extra large bore threaded connectors are female cap # 98713 and non-vented male cap. All connectors are available in custom colors for an additional charge.

Qosina has low minimums and products are in stock in our ISO 9001:2000 warehouse. Visit our website, [qosina.com](http://qosina.com) for a comprehensive list of over 5000 products. Call our Customer Specialists for assistance with your catalog/sampling request or order. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

