



**NEWS RELEASE  
FOR IMMEDIATE RELEASE**

## **Qosina Slide Clamps**

Edgewood, NY, March 28, 2005 – Qosina is proud to highlight our best-selling slide clamps. Qosina's slide clamps are composed of rigid high density polyethylene and fit 0.315 or 5/16 inch/8mm OD tubing. The clamps are 1-7/8 inches in length, 7/8 inch high and 1/8 inch thick. These clamps are available in several colors to help distinguish clamp locations.

- #11029 Red fits 0.315 or 5/16 inch/8mm OD tubing**
- #11031 Green fits 0.315 or 5/16 inch/8mm OD tubing**
- #11032 Light Blue fits 0.315 or 5/16 inch/8mm OD tubing**
- #11033 Blue Clamp fits 0.315 or 5/16 inch/8mm OD tubing**
- #11034 White Clamp fits 0.315 or 5/16 inch/8mm OD tubing**

Qosina has additional clamps, infusion and IV components to compliment your selections. We offer low minimums and products are in stock in our ISO 9001:2000 warehouse. Visit our website, [qosina.com](http://qosina.com) for a comprehensive list of over 5000 products. Call our Customer Specialists for assistance with your catalog/sampling request or order. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

