



NEWS RELEASE  
FOR IMMEDIATE RELEASE

## Qosina Sponge on a Stick

Edgewood, NY, May 17, 2005 –Qosina is pleased to add a new sponge on a stick, **Part # 13139** to our current inventory. The blue polyether foam head is approximately 1-1/2 inches by 1-1/2 inches by 7/8 inch thick. The ribbed polypropylene handle allows for a sure grip while the foam head, attached by adhesive, provides superior strength and durability. Overall length of this sponge on a stick is 6 inches.

Qosina offers many different sponges on applicator sticks in addition to various kit components. We offer low minimums and products are in stock in our ISO 9001:2000 registered warehouse. Visit our website, [qosina.com](http://qosina.com) for a comprehensive list of over 5000 products. Call our Customer Specialists for assistance with your catalog/sampling request or order. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

