

## Qosina Sterilization Indicators

Qosina has added more sterilization indicators to our ever-expanding product line. These indicators show a color change to signify proper exposure to various sterilization processes. Our indicators are available for:

- |  |        |                 |
|--|--------|-----------------|
| ◆ ETO ( <i>ethylene oxide</i> ) Sterilization Indicator        | #13122 | Purple to Green |
| ◆ Steam Sterilization Indicator                                | #13123 | Blue to Pink    |
| ◆ Radiation Sterilization Indicator ( <i>Gamma or E-Beam</i> ) | #13124 | Yellow to Red   |

These self-adhesive labels can be applied by hand. By the observable color change, it is easy to determine whether or not a particular sample has been through a qualified sterilization process. After exposure to the sterilization process, the indicators retain their color change indefinitely. The labels are approximately one-half inch in diameter and are available on rolls of 5,000.

Qosina offers additional sterilization accessories including CSR wraps and breathable pouches. Call our Customer Service Specialists who can assist you. To receive a sample, quotation or a copy of our 25<sup>th</sup> anniversary 440 page catalog, call (631) 242-3000, fax (631) 242-3230 or e-mail [info@qosina.com](mailto:info@qosina.com). You can view our extensive product line at [qosina.com](http://qosina.com). Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000 and ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

