

## Qosina Catheter Plugs

Edgewood, N.Y., June 3, 2005: Qosina is pleased to announce the addition of catheter plugs to our inventory. These plugs are ideal for quick and easy blocking of any size catheters. The tapered bodies push fit into tubing and have a flanged head to aid in removal. The polypropylene, dark blue plugs are available in three sizes:

		Length	Max. Diameter
◆ Small Catheter Plug	Part # 13131	56mm	7mm (21FR)
◆ Medium Catheter Plug	Part # 13132	67mm	9mm (27FR)
◆ Large Catheter Plug	Part # 13133	78mm	11mm (33FR)

To obtain a sample of our catheter plugs or other kit components call +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail [info@qosina.com](mailto:info@qosina.com). You can also request a free copy of our 25<sup>th</sup> Anniversary 440-page catalog with over 5000 SKUs. View our extensive product line at [qosina.com](http://qosina.com). Custom sourcing is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000 and ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

