



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

Qosina® Plastisol Reinforced Y Connector

Qosina is pleased to announce expanded production capabilities for one of our most popular components, **Part #81460 - Plastisol Reinforced Y Connector**. The connector is 2-3/8" (60 mm) long and has a bottom stem internal dimension of 1/4" (6.3 mm). The internal leg dimensions are undersized to insure a tight fit for OD tubing sizes of 3/32" (2.3 mm), 1/8" (3.1 mm), 3/16" (4.7 mm) and 1/4" (6.3 mm).

The plastisol connector is composed of USP Class VI liquid PVC resin. The beneficial feature of this universal connector is that it can be cut at various intervals to connect with assorted tubing sizes.

Qosina offers plastisol connectors in various sizes and configurations. We can create custom plastisol components to fit your requirements. Plastisol dipping can be implemented to create custom adapters, caps, connectors and other additional components. Custom multi-port Y connectors can be molded with as many as six ports to meet your design needs. Prototypes can be manufactured economically.

Please call our Customer Service Specialists who can assist you with your selection of connectors. To receive a sample, quotation or a copy of our 25th anniversary 440 page catalog, call (631) 242-3000, fax (631) 242-3230 or e-mail info@qosina.com. You can view our extensive product line at www.qosina.com. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

