



**News Release
For Immediate Release**

Qosina Tubing

Edgewood, N.Y., August 22, 2005 – Qosina now offers High Density Polyethylene (HDPE) tubing (**Part # T1150**) in coils of 250 feet (*76.2 meters*). The tubing has an ID of 0.10 inches (*2.54 millimeters*) and an OD of 0.15 inches (*3.81 millimeters*). This type of tubing is an excellent choice for storing and dispensing guidewires without kinking. Custom diameters are available and tubing can be cut to desired length.

Qosina has an extensive collection of tubing and guidewire accessories including j-straighteners, channel clips, torquers and bowls. Visit our website, qosina.com to view our complete offering of over 5,000 products. Call our Customer Specialists for assistance regarding a catalog, sampling request or an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.



QO-05-14-T1150