

Qosina Red Pilot Balloon with Plug

Edgewood, N.Y., October 4, 2005 – Qosina is pleased to add an additional pilot balloon to our current inventory. The red pilot balloon with plug - **part # 13499** is a component in respiratory and tracheal kits. The silicone balloon is 2-1/4 inches (57mm) in length with an 0.100 inch ID to accept the most popular inflation tubing. The port with the plug has an 0.158 inch ID. The balloons are brightly colored and can be easily distinguished from other components.

Additional products, tube components, inflation valves, luer activated checkvalves, connectors, suction ports and syringes can be purchased at Qosina. Call our Customer Service Specialists for assistance regarding obtaining a copy of our 2005 catalog, a sample request or to place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001:2004 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

