

Qosina Specimen Transport Bag

Edgewood, N.Y., October 4, 2005 – Qosina is adding a two-part low-density polyethylene (LDPE) specimen transport bag (Part # 13100) to our inventory. This specimen bag is leak resistant with two separate pockets usable for samples and documentation. The 6 inch by 9 inch (15.3cm by 23cm) bag features a fold-over Adhesive Ribbon™ seal. This creates a solid tamper-evident visual seal across both sections of the bag to keep employees and contents safe from contamination. These bags are available with custom absorbent inserts, printing and sizing options.

Qosina also offers additional sizes of specimen bags, transfer forceps, specimen cups and collection vials. Visit our website, qosina.com for a comprehensive listing of over 5,000 products. Call our Customer Specialists for assistance regarding a catalog, sampling request or an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is available by our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

