



**NEWS RELEASE  
FOR IMMEDIATE RELEASE**

## **Qosina Lighted Curettes**

Edgewood, N.Y., December 19, 2005 – Qosina is proud to present [single-use lighted curettes](#). Our unique curettes have a safe, internal battery powered LED light source providing full illumination in difficult to see locations. This convenient component permits single-handed use by eliminating the need for a separate electric light source. The light source attaches to the curette with a simple push and twist. Single-use curettes reduce cross-contamination and re-sterilization costs. The polycarbonate curettes are available with a choice of tips for a variety of applications.

- Q2210 [Lighted AngleLoop™ Curette](#), 50 pcs with 1 Light Source
- Q2220 [Lighted MicroLoop™ Curette](#), 50 pcs with 1 Light Source
- Q2230 [Lighted WaveCurette™](#), 50 pcs with 1 Light Source
- Q2240 [Lighted FlexLoop™ Curette](#), 50 pcs with 1 Light Source
- Q2260 [Lighted InfantLoop™ Curette](#), 50 pcs with 1 Light Source

Additional kit and procedure tray components are available from Qosina on our website at [www.qosina.com](http://www.qosina.com). Call our Customer Service Specialists to obtain a copy of our 2005 catalog, request a sample or to place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail [info@qosina.com](mailto:info@qosina.com). Custom sourcing service is available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001:2004 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

