



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

Qosina Multilayer Tubing

Edgewood, NY, July 20, 2006 — Qosina is pleased to add multilayer tubing ([**Part # T4200-4208**](#)) to our extensive inventory. This 78A durometer hybrid tubing combines the attributes of non-reactive polyethylene with clear, flexible medical grade PVC. The pure LDPE inner layer is compatible for contact with many reagents. The exterior non-DEHP PVC layer provides easy assembly, feel and bonding of traditional medical tubing. This tubing is available in coils of varying inner and outer diameters to meet your design needs.

| Part # | ID <i>in inches</i> | OD <i>in inches</i> |
|---------------|--------------------------------|--------------------------------|
| T4200 | .020" | .060" |
| T4201 | .042" | .090" |
| T4202 | .050" | .110" |
| T4203 | .060" | .125" |
| T4204 | .096" | .156" |
| T4205 | .120" | .160" |
| T4206 | .125" | .218" |
| T4207 | .156" | .218" |
| T4208 | .125" | .250" |

Qosina.com stocks many products to compliment your lines including luers, clamps, additional types of tubing and valves. Call our Customer Service Specialists to receive a copy of our new 2007 catalog, receive samples or to place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is also available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

