



**NEWS RELEASE  
FOR IMMEDIATE RELEASE**

## **Qosina Checkvalve Female Luer Inlet/Male Luer Outlet**

Edgewood, NY, May 26, 2006 –Qosina has added a new disc checkvalve ([Part # 11582](#)) with a female luer lock inlet and a male luer lock outlet. Our biocompatible checkvalve is comprised of lipid resistant, latex-free polycarbonate with a USP class VI silicone valve. This part is compatible with many compounds and has a 0.217 ppi cracking pressure. Our checkvalve features a color change from blue to clear occurring when sterilized using gamma radiation. This Qosina part is available with low minimums, providing an economical component for your needs.

Our checkvalves, tubing, spikes, luers, ports, bags and clamps in addition to tubing and other components are available to accommodate all your projects. Visit [qosina.com](http://qosina.com) or call our Customer Service Specialists to obtain our new 2007 catalog, request a sample or place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail [info@qosina.com](mailto:info@qosina.com). Custom sourcing service is available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001:2004 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

